



IR detection module

# **LAB Detector**

user's guide

# Table of Contents

1. Warranty.....	3
2. LAB Detector technical information.....	4
2.1. LAB Detector overview.....	4
2.2. LAB Detector parameters.....	4
2.3. Module physical dimensions.....	5
2.4. Absolute maximum ratings.....	6
2.5. Preamplifier load resistance.....	6
2.6. Supply socket.....	6
3. Getting started (regards the module itself as well as the complete IR detecting kit).....	7
3.1. Unpacking.....	7
3.2. Installation & connections.....	7
3.3. Switching on procedure.....	8
4. LAB Detector other related information.....	9
4.1. CE statement of conformity .....	9
4.2. Electromagnetic interference.....	9
4.3. Storage information.....	9
5. Safety instructions.....	10
6. Troubleshooting.....	11
6.1. List of the possible faults and solutions.....	11
6.1.1. Too high DC level at the output of the preamplifier.....	11
6.1.2. Too high current consumption by the preamplifier.....	11
6.1.3. No or too weak response to optical radiation.....	12
6.1.5. Excess noise.....	12
6.1.6. Unstable signal.....	12
6.1.7. Long signal rise/fall time.....	13

# 1. Warranty

VIGO System S.A. hereby represents and warrants all Products manufactured by VIGO and sold hereunder to be free from defects in workmanship or material during a period of twelve (12) months from the date of delivery save for products for which a special warranty is given. If any Product proves however to be defective in workmanship or material within the period herein provided VIGO System undertakes to the exclusion of any other remedy to repair or at its own option replace the defective Product or part thereof free of charge and otherwise on the same conditions as for the original Product or part without extension to original warranty time. Defective parts replaced in accordance with this clause shall be placed at the disposal of VIGO.

VIGO also warrants the quality of all repair and service works performed by its employees to products sold by it. In case the repair or service works should appear inadequate or faulty and should this cause malfunction or nonfunctioning of the product to which the service was performed VIGO shall at its free option either repair or have repaired or replace the product in question. The working hours used by employees of VIGO for such repair or replacement shall be free of charge to the client. This service warranty shall be valid for a period of six (6) months from the date the service measures were completed.

## **This warranty is however subject to following conditions:**

1. A substantiated written claim as to any alleged defects shall have been received by VIGO System within thirty (30) days after the defect or fault became known or occurred, and
2. The allegedly defective Product or part shall, should VIGO so require, be sent to the works of VIGO or to such other place as VIGO may indicate in writing, freight and insurance prepaid, properly packed and labeled.

## **This warranty does not however apply when the defect has been caused through**

1. normal wear and tear or accident;
2. misuse or other unsuitable or unauthorized use of the Product or negligence or error in storing, maintaining or in handling the Product or any equipment thereof;
3. wrong installation, assembly or failure to service the Product or otherwise follow VIGO's service instructions including any repairs or installation or assembly or service made by unauthorized personnel not approved by VIGO or replacements with parts not manufactured or supplied by VIGO;
4. modifications or changes of the Product as well as any adding to it without VIGO's prior authorization;
5. burned active element by irradiation above damage thresholds
6. electrostatic discharges
7. improper detector bias
8. improper TE cooler bias (TE cooler damage or active element overheating).
9. other factors dependent on the Customer or a third party.

Notwithstanding the aforesaid VIGO System liability under this clause shall not apply to any defects arising out of materials, designs or instructions provided by the Customer.

This warranty is expressly in lieu of and excludes all other conditions, warranties and liabilities, expressed or implied, whether under law, statute or otherwise, including without limitation any implied warranties of merchantability or fitness for a particular purpose and all other obligations and liabilities of VIGO or its representatives with respect to any defect or deficiency applicable to or resulting directly or indirectly from the Products supplied hereunder, which obligations and liabilities are hereby expressly canceled and waived. VIGO's liability shall under no circumstances exceed the invoice price of any Product for which a warranty claim is made, nor shall VIGO in any circumstances be liable for lost profits or other consequential loss whether direct or indirect or for special damage.

## **RMA Request Instructions:**

No Product may be returned without first contacting VIGO for a Return Material Authorization ('RMA') number. Please obtain a RMA number at [claim@vigo.com.pl](mailto:claim@vigo.com.pl) before returning any item. When requesting a RMA number, please state your order number, the product you wish to return and the reason for return. We will only accept returns which have a RMA number. Authorized returns are to be shipped according to received instruction from VIGO in appropriate shipping box. An unauthorized return, i.e. one for which an RMA number has not been issued and authorized returns however, shipped with incorrect customs documents - will not be accepted.

Please print the RMA number clearly on the return label to avoid any delay in processing. Please send package to:

**VIGO System S.A.**  
129/133 Poznanska St.,  
PL 05-850 Ozarow Mazowiecki  
Poland

## 2. LAB Detector technical information

Complete **LAB Detector** kit consists of:

- PVM-10.6 non cooled detector,
- VPDC-10I transimpedance preamplifier,
- PPS-02-12 power supply

### 2.1. LAB Detector overview

LAB Detector is a IR detecting kit delivered with a power supply and optional mounting post. Detector is sensitive to the radiation in the range from 2  $\mu\text{m}$  to 12  $\mu\text{m}$ . It makes it universal and sufficient to use in large variety of applications. Each preamplifier is also optimized for particular detector to achieve the best performance.

#### Applications:

- detection of cw, modulated and pulsed radiation within a broad 2 to 12  $\mu\text{m}$  spectral range
- pulsed, modulated & cw laser beam detection & monitoring
- laser beam diagnostics
- laser metrology
- laser - matter interaction studies
- laser based communication
- heterodyne detection


### 2.2. LAB Detector parameters

DETECTOR PVM-10.6			
Parameter	Unit	Value	Conditions
Spectral range	$\mu\text{m}$	2 to 12	
Time constant	nsec	$\leq 1$	
Detectivity	$\text{cm}\sqrt{\text{Hz/W}}$	$\geq 1 \cdot 10^7$	10.6 $\mu\text{m}$
Responsivity	V/W	$\geq 0.03$	10.6 $\mu\text{m}$
Detector area	$\text{mm}^2$	1 x 1	
Field of view	deg	90	

PREAMPLIFIER VPDC-10I			
Parameter	Unit	Value	Conditions
Bandwidth	MHz	10 MHz	-3 dB
Max. output signal amplitude	V	2	
Transimpedance	V/A	$1 \cdot 10^4$	$R_L = 1 \text{ M}\Omega$
Input current noise	pA/ $\sqrt{\text{Hz}}$	10	$F_o = 100 \text{ kHz}$
Output impedance	$\Omega$	50	
Detecting head diameter	mm	56	
Detecting head length	mm	108	
Weight	g	210	

PPS-02-12			
Parameter	Unit	Value	Conditions
Mains supply voltage	V	100 to 240	
Output voltage	V	+/- 12 V	$T_A = 22 \text{ }^\circ\text{C}$
Output current	mA	50	$T_A = 22 \text{ }^\circ\text{C}$

### 2.3. Module physical dimensions



Dimensions are given in millimeters.

## 2.4. Absolute maximum ratings

Parameter	immersed
Maximum Continuous Optical Power	2.5 W/cm <sup>2</sup>
Maximum Optical Power of Pulses Shorter Than 1µs	10 kW/cm <sup>2</sup>
Storage temperature range	-40 °C to +80 °C
Operating temperature range	-20 °C to +60 °C


## 2.5. Preamplifier load resistance

Although the output impedance is 50 Ω, preamp is suited to work with high impedance load.

**Attaching load with too low impedance may cause:**

- decreasing the transimpedance of the preamplifier,
- preamplifier overheat, especially if the output signal level is high,
- supply current abnormal increase.

## 2.6. Supply socket



**DB-9 connector**

Description	Pin number
GND	3
-V <sub>EE</sub>	6
+V <sub>CC</sub>	9

### 3. Getting started (regards the module itself as well as the complete IR detecting kit)

#### 3.1. Unpacking

Before unpacking, verify if there are no visible damages of the package. Report them to the manufacturer or the distributor.

#### 3.2. Installation & connections

- Mount the IR module in the measurement setup,
- Connect all elements with cables provided in the package,
- Make sure that all the elements are properly connected.

The picture shown below presents sample IR detecting kit configuration with interconnections.



### 3.3. Switching on procedure

- Do not switch on the preamplifier power supply until ALL elements of the IR detection kit are properly connected with the proper cables.
- Before starting the measurements make sure the load resistance is valid.
- Do not manipulate any connectors on the active IR kit – this may introduce additional noises to the system or can cause severe damages due to electrostatic discharges.



## 4. LAB Detector other related information

### 4.1. CE statement of conformity

VIGO System declares that IR detection kit is fully compliant with:

EU Low Voltage Directives: 73/23/EEC, 93/68/EEC

EU Electromagnetic Compatibility Directives: 89/336/EEC, 91/263/EEC, 92/31/EEC, 93/68/EEC

### 4.2. Electromagnetic interference

To avoid the electromagnetic interferences (EMI):

Please do not put any electromagnetic source close to the IR module. That may cause interference to the input of the transimpedance preamplifier and disturb the useful signal in the IR track















### 4.3. Storage information

To ensure reliable operation and long lifetime of detector, the following conditions should be fulfilled:

- Storage temperature:  $-10^{\circ}\text{C} \div +50^{\circ}\text{C}$  and 10% to 90% humidity,
- Avoid exposing to direct sunlight and strong UV/VIS light as this may result in degradation of a detector performance,
- Avoid electrostatic discharges (ESD),

## 5. Safety instructions

To ensure safe and failure-free operation of the IR detection kit comply with the following precautions:

-  The IR detector is dedicated to measure IR flux. Never use the detector in other applications.
-  Never dismount the case of the module. Do not open or remove the housing of the device. It may lead to its permanent damage or/and to a risk of electric shock.
-  Handle with care, avoid excessive mechanical shock or vibration. Excessive mechanical stress applied to the package itself or to a device containing the package may result in permanent damage.
-  Use the proper cables to connect elements of the kit, dedicated for specified device only. Never cut or shorten any cable, it may cause damage to your device.
-  Never use cables not approved by VIGO System.
-  Do not replace parts of the IR detection kit - cables, preamplifiers or detectors. Using parts from other IR detection kit can be the cause of malfunction of your module.
-  Never touch the detector window.
-  Do not insert any objects inwards the vents of the housing. Covering or closing the vents or cooling fan may cause a dangerous rise of the temperature inside the module.
-  All connectors must be connected carefully and firmly.
-  Avoid damp locations. The IR module should be neither installed nor used in high humidity locations.
-  Avoid locations with ambient temperature exceeding module operating temperature. The IR module should not be installed in the proximity of radiators, heat outlets, burners or other heat emitting objects.
-  Use proper ESD protection.
-  Ensure stable placement of the device. The IR module should stand on a stable surface.
-  Some IR detector window materials such as BaF<sub>2</sub> are soft. Particular attention should be paid not to scratch the surface of the window. A damaged window may entirely degrade the detector performance.

### Important note!

If you are setting up own system based on IR detecting system without power supply, please check if the preamplifier power supply are correct.

## 6. Troubleshooting

At first, please carefully read the previous parts of this manual, especially Safety instructions. If a problem persists, please check the information given below.

### General Test Procedure

#### Check:

- Whether components of your IR detection kit are connected correctly,
- Whether connections are snug-tight,
- Whether supply voltages or currents of the components are consistent with the specification,
- Attenuation/gain factors in your measurement set-up,
- Load resistance at the preamplifier output,
- For high frequency output line and the load must be matched to the 50 Ohm preamplifier output,
- Whether bandwidths of the components are reasonably well-matched to each other,
- Whether the detector is not saturated or damaged by too high irradiation. Refer to Beam Power Limitations section of this manual. If a problem still exists, try to identify and solve it by actions mentioned below,
- Cycle the power of the components off and on again,
- After shielding/switching off the radiation source, do you see any signal decrease? If yes, most probably you see a signal from the optical radiation you hope to see,

#### 6.1. List of the possible faults and solutions

##### 6.1.1. Too high DC level at the output of the preamplifier

Follow through General Test Procedure.

Remove the IR source from the detector field of view or decrease incident beam power of your IR source. If device still has high DC level, contact your technical support.

##### 6.1.2. Too high current consumption by the preamplifier

The device is probably damaged, contact your technical support.

### 6.1.3. No or too weak response to optical radiation

Follow through General Test Procedure.

- Check your radiation source,
- If you do not have any other detector to test your radiation source, move a hot object (ex. soldering iron) rapidly in front of the detector, across its whole field of view. (The higher the temperature of the object or shorter the distance to the detector, the higher the optical signal at the detector output.),
- Test your detector with another available radiation source. A warm or hot object moved in front of the detector can be used for DC or low frequency AC coupled devices,
- Align or improve your optical setup for the maximized signal or signal-to-noise ratio,
- Estimate the detector response, taking into account the incident radiation at the detector active area and the responsivity of the detector (or of the IR detector module), and check it against measured values,
- Evaluate the noise in the system. Calculate the noise introduced by the detector or detection module basing on their noise density from their data sheets. Check it against measured response of the blinded detector,
- Check if the incident radiation will be strong enough to obtain sufficient signal-to-noise ratio. (A signal is easy to pick out of a noise when its peak to peak value is several times greater than of the noise.) If it is lower than expected, consider increasing the incident radiation power on the detector active area,
- Check also if the incident radiation is not too strong, causing detector saturation or damage. If you cannot obtain signal sufficiently higher than noise, contact your technical support.

### 6.1.5. Excess noise

First of all check through General Test Procedure.

Excess noise may be caused by poor connections, ground connections and ground loops, high background photon flux or EMI (ex. inductive motors driving radiation chopper). You may also decrease the noise by reducing your system bandwidth.

If you are unable to identify external excess noise source, contact your technical support.

### 6.1.6. Unstable signal

Check through General Test Procedure.

Check:

- Cables,
- Radiation sources ,
- If the device is still unstable, contact your technical support.

Warning - DC or low frequency signal may vary due to fluctuations of thermal background radiation

### 6.1.7. Long signal rise/fall time

Check through General Test Procedure.

Check:

- The shape of the pulse with another detector if possible,
- The operating conditions (bias, optical setup),
- If impedance of the detector, cables, preamplifier and your read-out instrument are matched.